
BD

I attack
invaders **outside**
the cells.

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

B CELL


Summary of the key developmental stages and markers of B cells

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
E CD127 expression on Pro-B/Pre-B, Immature Transitional, Early Mature, and Late Mature B cells

Analysis of B-cell developmental stages in mouse bone marrow

C57BL/6 mouse bone marrow cells were stained with the following fluorescent antibodies: IgM FITC, CD43 APC, BP-1 PE, CD127 PE-Cy7, CD45R/B220 APC-Cy7, CD24 BD Horizon™ Brilliant Violet™ 421 (BV421), and IgD BD Horizon™ Brilliant Violet™ 510 (BV510), and analyzed using a BD FACSCanto™ II flow cytometer.

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
Backbone Markers	Human	Mouse
	CD19	CD45R/B220

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
Six-color analysis of human peripheral blood cells
Human peripheral blood mononuclear cells (PBMCs) were stained with the following fluorescent antibodies : CD19 PerCP-Cy[™]5.5, CD20 APC-H7, IgD FITC, CD27 PE-Cy[™]7, CD38 APC, and CD24 PE, and analyzed using a BD FACSVers[™] flow cytometer.

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Ten-color analysis of human peripheral blood cells
Human PBMCs were stained with the following fluorescent antibodies: CD19 BD Horizon[™] Brilliant Violet[™] 711, CD20 Alexa Fluor[®] 700, IgD PE-Cy7, CD27 BV421, CD38 APC, CD24 BD Horizon[™] PE-CF594, CD3 BD Horizon[™] V500, CD138 PE, IgM BD Horizon[™] Brilliant Violet[™] 605 (BV605), and IgG FITC, and analyzed using a BD LSRFortessa[™] flow cytometer.

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


Quantitation of immunoglobulins in human serum samples

Serum from five donors was tested using BD™ CBA Human Immunoglobulin Flex Sets to analyze total IgG, IgM, and IgA. Data was acquired on a BD FACSArray™ flow cytometer and analyzed using FCAP Array software. Dot plots for a negative control (detector alone) and one of the samples are shown. The results calculated (in g/L) on the basis of standard curves are shown, as well as the expected range.

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
Detection of Transcription Factors

Phospho-mTOR staining of human cells

Human peripheral B lymphocytes were stimulated with CpG ODN2395, then fixed with BD Cytotfix™ fixation buffer, permeabilized using BD Phosflow™ perm buffer III, and stained with anti-mTOR (pS2448) PE or a matching isotype control, and anti-CD20 Alexa Fluor® 647. Cells were analyzed using a BD FACSCanto™ II system.


Detection of Phosphoproteins

Blimp-1 staining in activated mouse splenocytes

B6 mouse splenocytes activated with LPS for 3 days were analyzed using the BD Pharmingen transcription factor buffer set, anti-CD45R/B220 PE, and either anti-Blimp-1 Alexa Fluor® 647 or a matching isotype control. Flow cytometry was performed using a BD FACSCanto II system.

Detection of Cytokines


XBP-1s expression in CpG-stimulated human PBMCs

CpG-stimulated human PBMCs were incubated with BD Horizon™ fixable viability stain 450, fixed and permeabilized using the BD Pharmingen transcription factor buffer set, and stained with anti-CD20 PerCP-Cy™ 5.5, and either anti-XBP-1s PE or a matching isotype control. Flow cytometry was performed using a BD FACSCanto II system.

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Comparison of technologies for T cell research

Tool/Technology	Flow Cytometry/Surface	Flow Cytometry/Intracellular	BD Cytometric Bead Array (CBA)	ELISPOT	ELISA	In Vivo Capture Assay
Molecules detected	Surface	Intracellular and surface	Secreted or intracellular	Secreted (in situ)	Secreted	Secreted (in vivo)
Multiparameter	Yes	Yes	Yes	No	No	No
Single cell/cell subset information	Yes	Yes	No	Frequencies, no subset information	No	No
Antigen specific	Yes	Yes	Yes	Yes	Yes	Yes
Post-assay viability	Yes	No	Yes, for secreted molecules	No	Yes	Yes
Quantitation of protein	Possible*	Possible*	Yes	No	Yes	Yes
Instrumentation	Flow cytometer	Flow cytometer	Flow cytometer	ELISPOT reader	Spectrophotometer	Spectrophotometer

*With a standard such as BD Quantibrite™ beads

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Major Known T-Cell Markers

Type of Cell	Cytotoxic	Th1	Th2	Th9 ¹	Th17	Tfh ²	Treg
Main Function	Kill virus-infected cells	Activate microbicidal function of infected macrophages, and help B cells to produce antibody	Help B cells and switch antibody isotype production	T cell proliferation and enhanced IgG and IgE production by B cells	Enhance neutrophil response	Regulate development of antigen specific B cell development and antibody production	Immune regulation
Pathogens Targeted	Viruses and some intracellular bacteria	Intracellular pathogens	Parasites	Parasites	Fungi and extracellular bacteria		
Harmful Function	Transplant rejection	Autoimmune disease	Allergy, asthma	Allergy	Organ-specific autoimmune disease	Autoimmune disease	Autoimmune disease, cancer
Extracellular Markers	CD8	CD4 CXCR3	CD4 CCR4, Ccrh2 (human)	CD4	CD4, CCR6	CD4, CXCR5	CD4, CD25
Differentiation Cytokines		IFN-γ, IL-2, IL-12, IL-18, IL-27	IL-4, IL-2, IL-33	IL-4, TGF-β	TGF-β, IL-6, IL-1, IL-21, IL-23	IL-12, IL-6	TGF-β, IL-12
Effector Cytokines	IFN-γ, TNF, LTα	IFN-γ, LTα, TNF	IL-4, IL-5, IL-6, IL-13	IL-9, IL-10	IL-17A, IL-17F, IL-21, IL-22, IL-26, TNF, CCL20	IL-21	TGF-β, IL-10
Transcription Factors		T-bet, Stat1, Stat6	GATA3, Stat6, Stat6	GATA3, Smads, Stat6	RORγt, RORγ, Stat3	Bcl-6, MAF	FoxP3, Smad3, Stat6


Markers can be altered as a result of cellular environment, differentiation state, and other factors. Key cytokines appear in bold. BD Biosciences offers reagents for molecules that are red.

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Marker	Fluorochrome	Purpose
Viability dye	V500	Viability
CD3	PerCP-Cy5.5	T cell marker
CD4	BUV395	T cell subsetting
CD8	FTC	T cell subsetting
CD127	Alexa Fluor® 647	Regulatory T Cell Marker
HLA-DR	APC-H7	Activation
CD45RO	PE-Cy7	Memory
CD197 (CCR7)	BV421	Naïve/Memory
CD38	BV605	Activation
CD27	BV786	Memory
CD25	PE-CF594	Regulatory T Cell Marker/Activation
CD196 (CCR6)	PE	Th17 Cell Marker
CXCR3	Alexa Fluor® 700	Th1 Cell Marker

Phenotyping of Cells with Unique Surface Profiles

The 13-color panel below examines memory and activation status of multiple T-cell subsets

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CAPABILITY	CBA	ICS
Allows detection of multiple cytokines in same experiment	✓	✓
Can obtain phenotype of specific cells expressing cytokine of interest		✓
Can measure quantity of cytokine secreted	✓	

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
Combining CBA and Intracellular Flow Cytometry to Examine Th17-cell Differentiation

With both CBA and intracellular cytokine staining (ICS) available, scientists at BD performed an experiment to examine T-cell differentiation, which can be induced by activation and treatment with cytokines. To study Th1/ Th2/Th17 cell differentiation, CD4⁺-panned human T cells isolated from normal donors were co-stimulated with CD3/CD28 and:


- IL-2, IL-12, and a neutralizing mAb to IL-4 (Th1 polarization)
- IL-2, IL-4, and a neutralizing mAb to IFN- γ (Th2 polarization)
- IL-2, IL-6, IL-1 β , TGF- β , IL-23, and a neutralizing mAb to IL-4 and IFN- γ (also tested with and without IL-2, IL-6, and TGF- β) (Th17 polarization)

Time Course Study of Th17 Cell Population

Representative data from Th17 polarized cell ICS experiments comparing levels of IL-17A with CD4, IFN- γ , and IL-4.

Cells were treated under the Th17 polarizing conditions described for the indicated time points. They were treated with BD GolgiStop (monensin) inhibitor, fixed and permeabilized with BD Cytofix/Cytoperm buffer, and then stained with antibodies against the indicated cytokines. At 6 days there were significant numbers of cells expressing IL-17A, with numbers of cells increasing at day 10 and then leveling off.

Data comparing cytokine levels as a result of different polarization conditions.


Supernatants from cells were polarized toward a Th1, Th2, or Th17 phenotype and cytokine levels were measured by CBA. As anticipated, each polarized condition resulted in the production of the signature cytokine associated with each Th cell type.

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Signaling pathway for Th17 cell differentiation

IL-6 stimulation leads to the phosphorylation at STAT3. STAT3 in concert with ROR γ T induced IL-17, the signature cytokine of Th17 cells.

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
T-cell activation profiles monitored with the BD Phosflow T Cell Activation Kit

The BD Phosflow™ Human T Cell Activation Kit is a comprehensive research system that uses flow cytometry to reliably determine the level of key phosphorylated signaling proteins involved in T-cell activation. The histogram overlays to the right show CD4⁺ and CD8⁺ T-cell signaling responses to treatment, monitored using the BD Phosflow Human T Cell Activation Kit. Response modifiers are described in the chart above.

BD

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