
Recent advances in single-cell sequencing hold great potential for exploring biological systems with unprecedented resolution.

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

How much RNA is in a single cell?

3 想做 Single Cell RNA-Seq, 但 Single Cell 內有多少 RNA 呢?

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
 Cryopreservation p.37





 Consumables p.38

isoplexis isoLight A Single-Cell Proteomics System

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IsoCode Chip

- Preinfusion polyfunctional anti-CD19 chimeric antigen receptor T cells are associated with clinical outcomes in NHL; Blood (2018) 132 (8): 804–814.
- Persistent Polyfunctional Chimeric Antigen Receptor T Cells That Target Glypican 3 Eliminate Orthotopic Hepatocellular Carcinomas in Mice; Gastroenterology, Volume 160, Issue 4, March 2021, Pages 1433
- TREATMENT USING CHIMERIC RECEPTOR T CELLS INCORPORATING OPTIMIZED POLYFUNCTIONAL T CELLS; WO/2018/187332


BD Rhapsody™ Single-Cell Analysis System

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
BD Rhapsody system workflow

BD Rhapsody Data

Severe COVID-19 Is Marked by a Dysregulated Myeloid Cell Compartment. Cell. 2020;182(6):1419-1440

Single-cell analysis reveals transcriptomic remodellings in distinct cell types that contribute to human prostate cancer progression. Nat Cell Biol. 2021;23(1):87-98.

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

- Cell count
- Cell viability

QC Report

- Cell/Bead loading efficiency
- Cell retention rate
- Bead retrieval efficiency
- Number of retrieved beads with viable cells

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

BD Rhapsody Publications

1. Modeling human adaptive immune responses with tonsil organoids. Nat Med. 2021;27(1):125-135.
2. Smc3 dosage regulates B cell transit through germinal centers and restricts their malignant transformation. Nat Immunol. 2021;22(2):240-253

portfolio


Quick selection guide

SMART-Seq[®] v4 Ultra[®] Low Input RNA Kit for Sequencing

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| Sequencing metrics comparing different cDNA synthesis protocols | | | | | | | |
|---|-----------|-----------------------------|--------|--------------|--------|------------|--------|
| RNA source | | 10 pg Mouse Brain Total RNA | | | | | |
| cDNA synthesis | | SMARTer ultra-low v3 | | SMART-Seq v4 | | SMART-Seq2 | |
| Yield (ng) | | 4.7 | 6.0 | 10.6 | 11.2 | 12.6 | 8.1 |
| Number of transcripts | FPKM >0.1 | 11,647 | 10,885 | 14,731 | 14,813 | 12,080 | 12,039 |
| | FPKM >1 | 9,729 | 9,105 | 12,501 | 12,591 | 10,270 | 10,058 |
| Percentage of reads (%): | | | | | | | |
| Mapped to rRNA | | 0.8 | 0.4 | 6.5 | 6.1 | 6.9 | 3.8 |
| Mapped to mitochondria | | 6.0 | 5.4 | 3.4 | 3.4 | 5.1 | 7.2 |
| Mapped to genome | | 96 | 97 | 96 | 95 | 72 | 93 |
| Mapped to exons | | 73 | 73 | 76 | 76 | 66 | 67 |
| Mapped to introns | | 21 | 21 | 19 | 20 | 28 | 27 |
| Mapped to intergenic regions | | 6.0 | 6.2 | 4.7 | 4.7 | 5.8 | 5.8 |

10 pg–10 ng of Human Brain Total RNA were used to generate cDNA libraries in duplicate with the SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing. cDNA libraries were amplified using 17, 14, 10, or 7 PCR cycles for the 10-pg, 100-pg, 1-ng, or 10-ng libraries, respectively. RNA-seq libraries were generated using the Nextera XT DNA Library Preparation Kit and sequenced on an Illumina MiSeq instrument. Sequences were analyzed as described in the Methods.

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
| Sequencing metrics comparing different cDNA synthesis protocols | | | | | | | |
|---|-----------|-----------------------------|--------|--------------|--------|------------|--------|
| RNA source | | 10 pg Mouse Brain Total RNA | | | | | |
| cDNA synthesis | | SMARTer ultra-low v3 | | SMART-Seq v4 | | SMART-Seq2 | |
| Yield (ng) | | 4.7 | 6.0 | 10.6 | 11.2 | 12.6 | 8.1 |
| Number of transcripts | FPKM >0.1 | 11,647 | 10,885 | 14,731 | 14,813 | 12,080 | 12,039 |
| | FPKM >1 | 9,729 | 9,105 | 12,501 | 12,591 | 10,270 | 10,058 |
| Percentage of reads (%): | | | | | | | |
| Mapped to rRNA | | 0.8 | 0.4 | 6.5 | 6.1 | 6.9 | 3.8 |
| Mapped to mitochondria | | 6.0 | 5.4 | 3.4 | 3.4 | 5.1 | 7.2 |
| Mapped to genome | | 96 | 97 | 96 | 95 | 72 | 93 |
| Mapped to exons | | 73 | 73 | 76 | 76 | 66 | 67 |
| Mapped to introns | | 21 | 21 | 19 | 20 | 28 | 27 |
| Mapped to intergenic regions | | 6.0 | 6.2 | 4.7 | 4.7 | 5.8 | 5.8 |

Replicate libraries were generated from 10 pg Mouse Brain Total RNA using the SMART-Seq v4 kit, the SMARTer Ultra low v3 kit, or the SMART-Seq2 method. 18 PCR cycles were used to amplify cDNA libraries with the SMART-Seq2 method and SMARTer ultra-low v3 kit; however, only 17 PCR cycles were needed for the SMART-Seq v4 libraries. RNA-seq libraries were generated using Nextera® XT DNA Library Preparation Kit and sequenced on an Illumina MiSeq® instrument.

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
Gene body coverage shown is the average of two replicate libraries prepared from 10 pg Mouse Brain Total RNA using the three different cDNA synthesis methods. The SMARTer Ultra low v3 kit produced a slight 3' bias, and the SMART-Seq v4 kit produced a slight 5' bias; however, the overall coverage was fairly even.

| 產品 | 包裝 | 貨號 |
|---|----------|--------|
| SMART-Seq® v4 Ultra® Low Input RNA Kit for Sequencing | 12 Rxns | 634888 |
| | 24 Rxns | 634889 |
| | 48 Rxns | 634890 |
| | 96 Rxns | 634891 |
| | 192 Rxns | 634892 |
| | 480 Rxns | 634893 |
| | 960 Rxns | 634894 |

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


Panel A. Libraries were prepared in duplicate from the same SSv4 cDNA using either the SMART-Seq Library Prep Kit portion of the SSv4 PLUS kit or Nextera XT. Panel B. Correlation plots between libraries made with SSv4 PLUS and Nextera XT kits. Pearson and Spearman's correlations were calculated from the TPM values obtained. Scatter plots showed the TPM values from all genes with a log₁₀+1 scale. Panel C. Venn diagram showing the number of genes identified from the processed cDNA, in duplicate, with the SSv4 PLUS kit, Nextera XT, or between both chemistries.

SMART-Seq v4 PLUS Kit 建構更高品質的 libraries

Ssv4 chemistry was used to produce cDNA, in triplicate, from 10 pg of Mouse Brain Total RNA. Illumina-compatible libraries were then generated from 1 ng or 125 pg using the SMART-Seq Library Prep Kit portion of the Ssv4 PLUS kit or Nextera XT, respectively, and sequenced on a NextSeq 500. The reads were normalized to 4M paired-end reads and analyzed as described in the methods. Panel B. The library yield obtained with the Ssv4 PLUS kit is higher than the yield obtained when using Nextera XT. Panel C. The number of genes identified was also higher for the Ssv4 PLUS kit compared to Nextera XT. Panel D. As expected, the distribution of the mapped reads is similar between the two library preparation methods. Panel E. Pearson and Spearman's correlations were calculated from the TPM values obtained from the triplicate cDNAs (A to C) processed with the Ssv4 PLUS or Nextera XT kit.

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

| Sequencing metrics comparing SMART-Seq v4 and SMART-Seq HT kits | | | | | | |
|---|-----------------------------|--------|--------|--------------|--------|--------|
| RNA source | 10 pg Mouse Brain Total RNA | | | | | |
| cDNA synthesis | SMART-Seq v4 | | | SMART-Seq HT | | |
| Replicate | A | B | C | A | B | C |
| Yield (ng) | 7.3 | 8.5 | 9.4 | 9.5 | 9.9 | 9.2 |
| Number of transcripts | 14,168 | 13,934 | 14,147 | 14,510 | 14,650 | 14,553 |
| | 11,455 | 11,267 | 11,349 | 11,582 | 11,670 | 11,454 |
| Average Pearson/Spearman | 0.97/0.67 | | | 0.97/0.68 | | |
| | 0.96/0.66 | | | | | |
| Proportion of reads mapped (%) | | | | | | |
| rRNA | 0.7 | 0.6 | 0.5 | 1.1 | 1.1 | 1.1 |
| Mitochondria | 2.8 | 4.2 | 4.1 | 3.7 | 3.8 | 4.2 |
| Genome | 92.3 | 90.8 | 86.6 | 89.1 | 89.0 | 87.8 |
| Exons | 74.8 | 72.5 | 69.0 | 67.3 | 67.1 | 65.3 |
| Introns | 13.3 | 14.1 | 13.3 | 16.9 | 16.9 | 17.4 |
| Intergenic regions | 4.2 | 4.2 | 4.4 | 4.9 | 5.0 | 5.2 |

Libraries were prepared from 10 pg of Mouse Brain Total RNA. The output cDNA was converted into RNA-seq libraries using the Illumina Nextera XT DNA Library Preparation Kit and sequenced on an Illumina NextSeq® instrument (2 x 75 bp). Sequences were analyzed as described in the methods after normalizing all the samples to 13

SMART-Seq HT PLUS Kit 建構更高品質的 libraries


SMART-Seq HT was used to produce cDNA from 10 pg of Mouse Brain Total RNA in triplicate. Illumina-compatible libraries were generated using SMART-Seq PLUS or Nextera XT (Illumina) library prep kits, and sequenced on a NextSeq 500 system. The reads were normalized to 4M paired-end reads and analyzed as described in the technical note.

| 產品 | 包裝 | 貨號 |
|-------------------------|----------|---------|
| SMART- Seq® HT Kit | 12 Rxns | 634455 |
| | 48 Rxns | 634456 |
| | 96 Rxns | 634437 |
| | 192 Rxns | 634438 |
| | 480 Rxns | 634436 |
| SMART- Seq® HT PLUS Kit | 48 Rxns | R400748 |
| | 96 Rxns | R400749 |

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

Single cells from the lymphoblastoid cell line GM12878 were processed with SSsc (18 cells) or the Smart-seq2 protocol (20 cells; Picelli et al. 2014) using 19 cycles of PCR. As described in the methods, RNA-seq libraries were generated and sequences analyzed (after normalizing all samples to 1.75 million paired-end reads).

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| Sequencing metrics comparing the SMART-Seq v4 kit and SMART-Seq Single Cell Kit | | | | | | |
|---|--------------------|-------|-------|-----------|--------|--------|
| RNA source | 2 pg UHR total RNA | | | | | |
| cDNA synthesis method | SSv4 | | | SSsc | | |
| Replicate | A | B | C | A | B | C |
| cDNA yield (ng) | 7.8 | 6.9 | 5.5 | 14.8 | 14.9 | 9.6 |
| Number of genes with TPM >1 | 7,412 | 7,522 | 7,487 | 8,774 | 8,614 | 8,406 |
| Number of genes with TPM >0.1 | 8,660 | 8,868 | 9,240 | 10,319 | 10,276 | 10,285 |
| Average Pearson/Spearman | 0.95/0.59 | | | 0.97/0.63 | | |
| Proportion of reads mapped (%): | | | | | | |
| Genome | 92.7 | 92.5 | 92.5 | 80.1 | 80.9 | 80.6 |
| Exon | 79.3 | 78.7 | 76.6 | 63.4 | 64.1 | 62.0 |
| Intron | 10.5 | 10.9 | 12.5 | 13.0 | 12.8 | 14.0 |
| Intergenic regions | 2.9 | 3.0 | 3.4 | 3.7 | 4.0 | 4.6 |
| rRNA | 0.8 | 0.7 | 0.6 | 6.1 | 6.0 | 4.3 |
| Mitochondria | 3.5 | 3.6 | 3.9 | 9.3 | 8.4 | 10.2 |


Replicate cDNA libraries were generated from 2 pg of Universal Human Reference (UHR) total RNA using the SMART-Seq v4 kit (SSv4) or the SMART-Seq Single Cell Kit (SSsc); all libraries were processed with 19 PCR cycles. As described in the methods, RNA-seq libraries were generated and sequences analyzed (after normalizing all samples to 1.6 million paired-end reads). SSsc identified about 15% more genes than SSv4

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
Libraries made from 2 pg of UHR total RNA (Table I) were analyzed using scatter plots to visualize the reproducibility between technical replicates (shown are TPM values from all genes with a log10+1 scale). SSV4 (Panel A) generated highly reproducible quantification, but SSsc (Panel B) produced superior reproducibility, as seen in the increased Pearson and Spearman correlations. In addition, SSsc was better at detecting low-expression genes.

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
12 single cells from lymphoblastoid cell line GM22601 were processed with SSV4 or SSsc using 19 cycles of PCR. As described in the methods, RNA-seq libraries were generated and sequences analyzed (after normalizing all samples to 1.25 million paired-end reads). For all boxplots in this figure, the box denotes the interquartile range (IQR), i.e., the 25th and 75th quartiles; the whiskers are 1.5 x IQR from the median value and represent the extremes of the data.

增進 primary samples 的定序表現

PBMCs from a healthy donor were processed with the SSv4 or SSsc kit (~50 single cells per kit). RNA-seq libraries were generated as described in the methods. Panel A. The read distribution is fairly similar between the two chemistries. Panel B. About 60% more genes are detected in the cells processed with SSsc, regardless of the number of reads used for the analysis.


Libraries were prepared from T cells according to manufacturer's instructions with either SSsc or NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina. Panel A. In these boxplots, the box denotes the interquartile range (IQR), i.e., the 25th and 75th quartiles; the whiskers are 1.5 x IQR from the median value and represent the extremes of the data. The number of genes identified with a TPM > 0.1 is higher (~40%) for SSsc (median = 3,591) than for NEBNext (median = 2,665). Panel B. The read distribution is different between the two chemistries, with more reads mapping to exon regions for SSsc. Panel C. The high number of cycles required to generate a library for the NEBNext protocol means that the negative controls for NEB (red) cluster strongly with the single-cell libraries (light blue), while the negative controls for SSsc (purple) are distinct from the single-cell libraries (dark blue).

| 產品 | 包裝 | 貨號 |
|----------------------------|---------|--------|
| SMART-Seq® Single Cell Kit | 12 Rxns | 634470 |
| | 48 Rxns | 634471 |
| | 96 Rxns | 634472 |

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



SMART-Seq Stranded Kit 彈性的樣本起始量

| Sequencing alignment metrics for A375 total RNA and cells | | | | | | | |
|---|-----------|-------------|-----------|-----------|-----------|-----------|-----------|
| Input | Total RNA | 1,000 cells | 500 cells | 100 cells | 10 cells | 5 cells | 1 cell |
| Number of reads (pairs) | 6,000,000 | 6,000,000 | 6,000,000 | 6,000,000 | 6,000,000 | 6,000,000 | 5,873,974 |
| Number of transcripts >1 FPKM | 13,260 | 13,294 | 13,583 | 13,520 | 12,726 | 12,602 | 11,540 |
| Number of transcripts >0.1 FPKM | 21,334 | 21,113 | 21,365 | 21,145 | 20,550 | 18,888 | 15,815 |
| Proportion of reads (%): | | | | | | | |
| <i>Exonic</i> | 34.7 | 36.4 | 39.2 | 42.7 | 36.7 | 36.2 | 37.3 |
| <i>Intronic</i> | 29.6 | 29.3 | 27.7 | 28.3 | 34.0 | 30.4 | 21.1 |
| <i>Intergenic</i> | 14.2 | 13.4 | 12.2 | 12.9 | 16.7 | 16.8 | 10.1 |
| <i>rRNA</i> | 7.0 | 11.4 | 11.5 | 6.3 | 3.6 | 4.9 | 7.1 |
| <i>Mitochondrial</i> | 4.1 | 3.5 | 3.7 | 4.9 | 3.8 | 4.4 | 4.6 |
| Overall mapping (%) | 89.6 | 93.9 | 94.3 | 95.1 | 94.9 | 92.7 | 80.2 |
| Duplicate rate (%) | 37.3 | 45.2 | 40.3 | 46.1 | 52.5 | 72.2 | 78.5 |
| lncRNA mapping: | | | | | | | |
| Number of mapped reads (%) | 7.2 | 10.4 | 10.8 | 9.4 | 8.7 | 8.6 | 7.3 |
| lncRNA transcripts detected | 5,395 | 4,687 | 4,565 | 5,439 | 5,440 | 4,983 | 2,802 |

A375 cells isolated by FACS were used to generate RNA-seq libraries with the SMART-Seq Stranded Kit. Input varied from 1 cell to 1,000 cells, with two replicates per input of 5–1,000 cells and 12 replicates for the single cells. For comparison, two aliquots of 1,000 cells were used for total RNA purification and then used for library preparation.

SMART-Seq Stranded Kit 與 SMART-Seq v4 Kit 具有一致性的靈敏度


Single cells (K562) isolated by FACS were used to generate RNA-seq libraries with the SMART-Seq Stranded Kit (Stranded) and a SMART-Seq v4 Kit (SSv4; cDNA from this kit was further processed with a Nextera® XT DNA Library Preparation Kit). Panel A. The overlap in the total number of transcripts identified (FPKM >1) by each kit was analyzed and shown to be 83.8%. Panel B. The number of transcripts identified (FPKM >1) in individual cells was similar between the two kits, with a tighter range across cells processed with the SSv4 kit. Panel C. The reproducibility (Pearson correlation) of transcript expression levels across all cells from each kit was similar, although slightly higher and more consistent across cells processed with the Stranded kit.

| 產品 | 包裝 | 貨號 |
|-------------------------|---------|--------|
| SMART-Seq® Stranded Kit | 12 Rxns | 634442 |
| | 48 Rxns | 634443 |
| | 96 Rxns | 634444 |


SMARTer[®] Human scTCR a/b Profiling Kit

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
Panel A. Schematic of a T-cell receptor comprised of an alpha chain (TCR- α) and a beta chain (TCR- β). Panel B. Schematic representing the difficulties of obtaining pairing information for alpha and beta chains from bulk sequencing data. Rare clonotypes can allow for assessing the pairing of some cells (clonotypes in orange).



Pooling 策略

SMART-Seq Indexed Oligo A
 SMART-Seq Indexed Oligo B
 SMART-Seq Indexed Oligo C
 SMART-Seq Indexed Oligo D
 SMART-Seq Indexed Oligo E
 SMART-Seq Indexed Oligo F
 SMART-Seq Indexed Oligo G
 SMART-Seq Indexed Oligo H

Samples are pooled by column, such that each pool contains eight cells, each with a differently indexed SMART-Seq Indexed Oligo. Different combinations of the Forward and Reverse HT indexes are used during PCR 2 to allow multiplexing of the samples into a single flow-cell lane (see the User Manual for more details).

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Panel A. Cell-type calling based on the identified clonotypes for each well. The seven omitted cells did not have clonotype calls for either TCR-α or TCR-β with read numbers that were above the threshold. Panel B. Analysis of pairing information. Paired TCR-αβ clonotypes were obtained for 34 cells in the plate. Panel C. The alpha, beta, or alpha-beta pairing information represented as a percent of cells analyzed. Omitted cells were not included in this analysis.


| 產品 | 包裝 | 貨號 |
|--|----------|--------|
| SMARTer® Human scTCR a/b Profiling Kit | 96 Rxns | 634431 |
| | 480 Rxns | 634432 |

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



Flow-sorted cancer cells were amplified with PicoPLEX single cell WGA technology. Each of the cells (blue lines) amplified at the same rate and resulted in a similar, predictable yield. A pooled sample of five cells (red line) represents a positive control. Samples containing 0 cells (green lines; no template control) show very low background.

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Single-blastomere biopsies were amplified using PicoPLEX single cell WGA technology, labeled, and hybridized to BlueGnome's 24Sure arrays at Genesis Genetics Institute. Note the clear indication of CNVs. In 2011, ESHRE clinical trials confirmed accuracy of karyotyping using PicoPLEX technology.

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Locus-specific qPCR was used to quantify 48 loci in independent single-cell libraries. Data shown compares results of DNA from two individual samples amplified with (left to right) PicoPLEX single cell WGA technology, GenomePLEX, and GenomiPhi using 10 pg of DNA. More than 90% of the product from the PicoPLEX single cell WGA technology produced identifiable, highly reproducible human sequences.

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
The PicoPLEX WGA v3 system is more accurate at detecting SNVs when compared to DOPlify and REPLI-g technologies

| Depth of SNV position ≥ 10 Allele frequency $\geq 20\%$ | PicoPLEX WGA v3 1 cell | | PicoPLEX WGA v3 5 cells | | DOPlify 1 cell | | DOPlify 5 cells | | REPLI-g 1 cell | | REPLI-g 5 cells | | |
|---|------------------------|-------|-------------------------|--------|----------------|-------|-----------------|--------|----------------|--------|-----------------|-------|--------|
| | Bulk | Rep1 | Rep2 | Rep1 | Rep2 | Rep1 | Rep2 | Rep1 | Rep2 | Rep1 | Rep2 | Rep1 | Rep2 |
| Number of SNVs called | 74 | 57 | 67 | 69 | 67 | 34 | 57 | 62 | 67 | Failed | Failed | 40 | Failed |
| Number of false positives | | 3 | 1 | 0 | 1 | 5 | 1 | 0 | 7 | Failed | Failed | 0 | Failed |
| Average false positives | | 0.02% | | 0.005% | | 0.03% | | 0.035% | | Failed | | 0 | Failed |
| Call rate | | 78% | 92% | 95% | 92% | 47% | 78% | 85% | 92% | Failed | Failed | 55% | Failed |
| Average call rate | | 85% | | 93% | | 62% | | 88% | | Failed | | 55% | Failed |
| Missed | | 17 | 7 | 5 | 7 | 40 | 17 | 12 | 7 | Failed | Failed | 34 | Failed |
| Average locus dropouts | | 16.2% | | 8.1% | | 38.5% | | 12.8% | | Failed | | 45.9% | Failed |
| Number of heterozygous SNVs called | 45 | 45 | 38 | 45 | 45 | 36 | 31 | 38 | 41 | Failed | Failed | 32 | Failed |
| Average allele dropouts | | 7.8% | | 0.0% | | 25.6% | | 12.2% | | Failed | | 71.1% | Failed |


Whole genome amplification products from single- or five-cell samples of a GM12878 cell line (Coriell Institute) were prepared in replicates, and SNVs detected were reported as numbers and percentages. Although the REPLI-g system produced sufficient yield, only one of the five-cell samples contained enough amplicon material to sequence, and, therefore, no data is available for the other samples. An intersection of Genome In a Bottle (GIB) variants to hg19 (human genome assembly GRCh37, Ensembl) indicates that a total of 78 variants are expected to be present for the GM12879 cell line. Due to the amplicon design, paired-read lengths of 75 bp were too short to capture 4 out of the 78 SNVs; therefore, the total number of capturable SNVs was cut down to 74. VarDict was used to interpret SNVs from BAMfiles using the following criteria: depth of SNV position ≥ 10 reads (10X coverage), allele frequency $\geq 20\%$.

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Single cells from various cell lines (GM22601, GM05067, and GM12878) were amplified using a prototype of the PicoPLEX Single Cell WGA Kit v3. 1 ng of the purified product was used as input for a Nextera® XT kit and sequenced on an Illumina MiSeq® platform at a read length of 2 x 75 bp. Fastq files were trimmed to remove adapters and then aligned to the human genome assembly GRCh37. Only autosomes are reported. For both panels, alignment was normalized to 1 million reads, and the number of reads per bin of 1 Mb was calculated using bedtools 2.25.0. The log2 ratio of the bin counts (Sample/Reference) was plotted using the Integrative Genomics Viewer.

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
Whole genome coverage reproducibility of the prototype PicoPLEX Single Cell WGA Kit v3 in comparison to DOPlify and REPLI-g (MDA) kits. WGA products were prepared from single-cell samples of GM12878 in replicates, using a prototype of PicoPLEXWGA v3, DOPlify, and REPLI-g (MDA) kits. 1 ng of amplified product was used as input for a Nextera XT kit, and the resulting libraries were sequenced on an Illumina MiSeq platform using a read length of 2 x 75 bp. After read alignment to human genome assembly GRCh37 and normalization to 1 million reads (0.5 million read pairs), the number of reads per bin of 1 Mb was calculated using bedtools 2.25.0. Total reads in each window from two single-cell libraries were plotted, with Pearson and Spearman correlations calculated and indicated on each graph

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


Panel A. A log-log plot showing the number of bases covered at various depths of sequencing (~35M read pairs, PE 2 x 150 bp). The coverage of PicoPLEX Gold was similar to QIAseq FX (MDA) at lower depths and greater at higher depths. Panel B. The reproducibility of coverage evaluated by comparing total reads in 100-kb bins. The consistency of the total reads in each window from the two single-cell libraries is significantly higher for PicoPLEX Gold (left), both in comparison to QIAseq FX (right) and to other technologies in the market (data not shown).

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Panel A. Comparison of SNV-detection rate between QIAseq FX (MDA), PicoPLEX Gold (PP Gold), and SMARTer PicoPLEX WGA (PP WGA) kits from one cell (1c), five cells (5c) or 15 pg of NA12878 gDNA inputs. Single and five cells were sequenced to a depth of ~35M read pairs, and gDNA samples to a depth of ~40M read pairs. The high fidelity and robust coverage of PicoPLEX Gold (blue bars) provide a clear advantage in detecting a greater (~2–9 fold) number of high-quality SNVs compared to QIAseq FX (gray bars) and PicoPLEX WGA (purple bar). Panel B. The symmetric distribution of the B-allele frequencies for PicoPLEX Gold (blue bars), centered around 0.5, indicating a balanced recovery of both alleles. PicoPLEX has better allele balance compared to QIAseq FX (MDA) (gray bars). Panel C. Unbiased amplification of PicoPLEX Gold results in the lowest allele drop-out (false-negative) rates among all single cell library-preparation technologies tested. Panel D. High fidelity of the polymerases used in PicoPLEX Gold kit (blue bars) leads to minimal allele drop-in rates that are comparable to QIAseq FX (gray bars) and significantly lower than PicoPLEX WGA (purple bar).

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Log₂ ratio of the total number of reads in 50-kb bins from single NCI-H929 cells, shown as one cell in Panel A and a second cell in Panel B. Red bars represent copy-number gains while blue bars represent losses. The top row of the graphs in each panel depicts the control bulk sample sequenced to a depth of 90 million read pairs. The highly reproducible coverage of the SMARTer PicoPLEX Gold kit enables the accurate detection of structural variants as small as 100 kb, even at shallow sequencing depths (2.5–8.5 million read pairs).

| 產品 | 包裝 | 貨號 |
|---|----------|---------|
| SMARTer® PicoPLEX® Gold Single Cell DNA -s eq Kit | 24 Rxns | R300669 |
| | 96 Rxns | R300670 |
| | 384 Rxns | R300698 |

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| Cell type | Description | | |
|---------------|--|----------|--|
| P3/x63-Ag8.U1 | Murine myeloma cell | SK007 | Human B-cell line |
| 2D-8 | Murine hybridoma | K562 | Human caucasian chronic myelogenous leukaemia cell |
| YAC-1 | Murine lymphoblast | HeLa | Human uterine cervical carcinoma cell |
| NBM-Lu | Normal newborn murine fibroblast cell line | HepG2 | Human hepatocellular carcinoma cell |
| 129SV | Mouse ES cell | Caco-2 | Human colonic adenocarcinoma cell |
| Feline PBMC | Feline peripheral blood mononuclear cells | UE6E7-16 | Human Mesenchymal cell |
| Canine PBMC | Canine peripheral blood mononuclear cells | UE7T-13 | Human Mesenchymal Stem cell |
| Jurkat | Human T-cell line | 201B7 | Human iPSc cell |

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Technology Inc.