immunohistochemistry

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0800-211-819

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

小鼠胸線組織以 neutral buffered formalin 固定

IHC Zinc Fixative: Mouse thymus fixed in IHC Zinc Fixative (left panel) or neutral buffered formalin (right panel) was embedded in paraffin and sectioned. The sections were stained with Purified Mouse anti-Rat Crry/p65 antibody (Cat. No.554991)

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

BD Retrievacen 可優化 IHC 染色結果

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

Antigen retrieval by BD Pharmingen Retrievacen A (pH 6.0): Formalin-fixed paraffin sections of normal human tonsil were either untreated (left panel) or pre-treated with BD Retrievacen A (pH 6.0) (right panel) and then stained with the anti-Ki-67 antibody. Note the absence of antibody staining in the untreated versus the panel where antigen retrieval solution was used.

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Antigen retrieval by BD Retrievacen B (pH 9.5): Formalin-fixed paraffin sections of normal human tonsil were either untreated (left panel) or pre-treated with BD Retrievacen B (pH 9.5) (right panel) and then stained with the anti-human topoisomerase II antibody. Note the absence of antibody staining in the untreated versus the panel where antigen retrieval solution was used.

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
Immunostaining with antibodies against MLH1 and PMS2. Thick arrows indicate tumor cells. Thin arrows indicate internal control cells, either stromal or epithelial. A and B, well-differentiated MSI-H colon carcinoma (hMLH1 germ-line mutation, exon 16 delK618) with negative MLH1 and PMS2 staining in tumor nuclei, with retained staining of stromal cells.

參考文獻：Clin Cancer Res. 2004 Feb 1;10(3):972-80. Microsatellite instability, immunohistochemistry, and additional PMS2 staining in suspected hereditary nonpolyposis colorectal cancer.

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
Characterization of cell lines and fragment-derived tumours from mouse angiosarcomas. The morphology of tumours developing following implantation of angiosarcoma-derived cell lines was that of undifferentiated pleomorphic sarcomas. H&E staining (i) and immunohistochemical analysis using antibody to CD31 (ii). The tumour cells were negative for CD31, as was the implanted angiosarcoma-derived cell line (inset).

參考文獻 : Dis Model Mech. 2019 Jul 9;12(7). Development of mouse models of angiosarcoma driven by p53.

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Mouse/Rat/Hamster Ig HRP Detection Kits

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- Biotinylated anti-Ig secondary antibody
- Antibody diluent
- Pre-diluted Streptavidin-HRP
- DAB substrate

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
- 適用於 acetone-fixed frozen sections, formalin- or zinc-fixed paraffin sections, cytopins/smears
- 通用款，適合任一種 IHC purified primary antibody
- 專一性高，降低背景值的干擾
- 提供 optimized protocol，方便，增加實驗再現性

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
Immunohistochemical staining for MMP-2. Paraffin sections of human breast cancer tissue were stained with the A-Gel VC2 antibody specific for human MMP-2. The BD Pharmingen™ Anti-Mouse Ig HRP Detection Kit (Cat. No. 551011) was used to develop the signal. Cells positive for MMP-2 can be identified by the brown staining. Magnification 20 $\times$.

Immunohistochemical staining of CD8+ T lymphocytes. Frozen sections of normal mouse spleen were stained with 53-6.7 antibody specific for mouse CD8a. The BD Pharmingen™ Anti-Rat Ig HRP Detection Kit (Cat. No. 551013) was used to develop the signal. CD8+ T lymphocytes can be identified by the intense brown labeling of their cell-surface membranes. Magnification 20 $\times$.

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抗體分類	貨號	相對應抗原	Clone	樣本物種	包裝
DNA repair	556415	PMS2	A16-4	Hu, Ms	0.1 mg
Cell adhesion	550274	CD31	MEC13.3	Ms	1 mL
Stem cell	550538	CD44	IM7	Ms	1 mL
	550392	CD44	G44-26	Hu	1 mL
	550402	CD90	5E10	Hu	1 mL
	550543	CD90.2	53-2.1	Ms	1 mL
	550546	CD105	MJ7/18	Ms	1 mL
	550390	CD34	581	Hu	1 mL
Immune	550374	CD11b/Mac-1	ICRF44	Hu	1 mL
	550375	CD11c	B-ly6	Hu	1 mL
	550376	CD14	M5E2	Hu	1 mL
	550368	CD3	UCHT1	Hu	1 mL
	550372	CD8	HIT8a	Hu	1 mL
	550282	CD11b	M1/70	Ms	1 mL
	550280	CD4	RM4-5	Ms	1 mL
	550281	CD8a	53-6.7	Ms	1 mL
	550550	H-2Kb	AF6-88.5	Ms	1 mL

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