




alvetex 3D cell culture

讓每個細胞生長宛如組織，
維持生理上的 3D 形狀和結構

TECHNOLOGY TO ENABLE 3D CELL CULTURE TO BE ACHIEVED SIMPLY AND ROUTINELY



Ordering Information

Product Name	Code	Description
alvetex® 12 well plate	AVP002	1 x alvetex® 12 well plate
alvetex® 6 well inserts	AVP004-3	3 x individually sealed 6 well inserts
alvetex® 12 well inserts	AVP005 3	3 x individually sealed 12 well inserts
alvetex® well insert holder and petri-dish	AVP015	1 x well insert holder in a deep petri-dish with lid

Synthetic retinoid ~ ec23

穩定且可控制 stem cell 及 progenitor cells 生長分化的小分子藥物。

- 1 為 pan-RAR receptor agonist，具有與 all-trans retinoid acid (ATRA) 相同的生物活性
- 2 相較於 all-trans retinoid acid，ec23 對於光及熱較不敏感，能在一般實驗室的環境下保有一定的物理及化學特性



Neurite outgrowth from neuroprogenitor cells treated with ec23® for 21 days