

獨家Premix設計 · 限量體驗



- Premix 形式：減少污染及試劑的浪費
- 靈敏度高
- 再現性高：減少 pipetting error，增加實驗再現性
- 適用性廣：PVDF 及 NC 皆適用
- 選擇性多：針對不同的實驗需求，提供不同靈敏度等級的呈色劑

靈敏度
高

再現性
高

適用性
廣

選擇性
多

超激

Luminata Classico

Luminata Crescendo

Luminata Forte



Unique Feature	Premixed	Premixed	Premixed
Best used for	Blots where the primary antibody is incubated for ~1 hr	Blot where the primary antibody is incubated > 2 hrs	Blots with overnight primary antibody incubation, or detection of PTM** proteins
Detection Range	~6 pg	~1-3 pg	~400 fg
Signal Duration	1 hr	3 hr	3 hr
Stock Solution Stability	1 yr at 4 °C	1 yr at 4 °C	1 yr at room temperature

**PTM – Post-translationally modified

超低特價
方案



分類	Cat No.	Description	Qty/Pk	超值價
低階	WBLUC0500	Luminata Classico Western HRP substrate	500 mL	3,950
中階	WBLUR0500	Luminata Crescendo Western HRP substrate	500 mL	5,000

超值加購
方案



凡原價購買 Merck & Millipore 抗體 2 支
即可享 Luminata Substrate 加購價

Luminata Substrate
加購價

5折